**REGENFAST®**

 **Mastelli**

confinis ch-rep ag  
Hauptstrasse 16  
3186 Düringen / Switzerland

Geistlich Pharma AG,  
Bahnhofstrasse 40, 6110  
Wolhusen, Switzerland.  
Tel.: +41 41 492 55 55

UK RP: Geistlich Sons Limited  
Thorley House,  
Bailey Lane, Manchester Airport,  
GB - Manchester M90 4AB

Mastelli srl  
Via Bussana Vecchia,32  
18038 Sanremo (IM) Italy  
+39 0184 5111  
[www.mastelli.it](http://www.mastelli.it)

CE 0373

# REGENFAST

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05/2023

# REGENFAST

Sterile, viscoelastic gel with polynucleotides (10 mg/mL) and hyaluronic acid (10 mg/mL) for use in dentistry

EN

## PRODUCT DESCRIPTION

REGENFAST is a sterile, viscoelastic gel in a disposable pre-filled syringe for use in dentistry. REGENFAST benefits the enhancement of the microenvironment and cell vitality and improves tissue trophism by promoting more efficient tissue repair. Bioregeneration is the strengthening and optimisation of the physiological activities of cells involved in the natural regenerative processes of tissues in the human body, thanks to the action of organic molecules such as polynucleotides.

The polynucleotides contained in the device are highly purified, completely resorbable and natural substances derived from fish. The hyaluronic acid contained in the device is not of animal origin but is bio-fermented. It is one of the fundamental components of the extracellular matrix and imparts mechanical and elastic properties to the tissues.

## INDICATIONS

REGENFAST with polynucleotides and hyaluronic acid is indicated for the protection of the tissues of the oral cavity and the promotion of faster, physiological healing. REGENFAST can be used in the oral environment for surgical applications.

REGENFAST is intended for applications in periodontal defects in the context of surgical and non-surgical procedures. It can be mixed with (autologous) bone or bone substitutes to foster their integration in the context of regenerative therapy.

REGENFAST is intended for local application to improve trophism of the gum tissues and interdental papilla on attached and free gingiva.

## COMPOSITION

Each pre-filled syringe vial contains isotonic gel comprising: polynucleotides (10 mg/mL) and hyaluronic acid sodium salt (10 mg/mL), mannitol, sodium phosphate monobasic dihydrate, sodium phosphate dibasic dodecahydrate, water for injection.

## INSTRUCTIONS FOR USE

REGENFAST is a medical device in a pre-filled syringe vial containing a ready-to-use sterile gel. Prior to treatment, the user must inform their patient of the indications, contraindications, interactions, potential side effects and risks associated with the use of the product and must ensure that the patient is aware of signs and symptoms that may indicate potential complications.

Before starting treatment, the periodontium must be pre-treated by causal or ultrasound therapy to remove purulent material or necrotic debris. REGENFAST can also be used for subgingival applications. REGENFAST should be spread evenly over the entire surface of the affected area, if necessary using the most suitable applicator (usually a blunt tip needle).

The pre-filled syringe vial is ready for use: simply remove the protective cap and screw on the most suitable applicator.

The dosage should be adapted to the patient's clinical picture.

## **CONTRAINDICATIONS**

Use of REGENFAST is contraindicated in patients with individual hypersensitivity to the product's ingredients or with a history of allergy to fish products. Because no clinical data are available on the use of the product by individuals with autoimmune diseases, by pregnant or lactating women or by persons under the age of eighteen, treatment with REGENFAST in these cases is not recommended.

## **PRECAUTIONS**

Before, after, and during treatment with REGENFAST, observe standard rules of dental hygiene and prophylaxis.

Thoroughly disinfect the area to be treated before applying the product. In the event of an infection, precede application of REGENFAST with a specific treatment.

## **WARNINGS**

REGENFAST should only be administered by a specialist doctor or dentist or by personnel authorised to do so in compliance with local legislation and who are familiar with the correct procedures for carrying out the various stages of treatment.

Keep out of reach of children.

Check the integrity of the packaging and syringe before use.

Do not use this product if the packaging has been damaged.

Do not use this product after the expiry date indicated on the packaging. The expiry date refers to correctly stored intact product.

Once opened, the product must be used immediately and any residual contents must be discarded. This is a disposable device. Do not re-sterilise. Dispose of the product responsibly after use. Discard the used product in compliance with medical and biological waste disposal regulations. Do not freeze.

## **INTERACTIONS**

As hyaluronic acid sodium salt can precipitate in the presence of quaternary ammonium salts, the simultaneous use of REGENFAST and any products containing such substances, is to be avoided. It is therefore recommended that before application any disinfectant be removed using sterile saline solution.

There are currently no known possible interactions of the device with other devices, medicines or substances.

## **STORAGE INSTRUCTIONS**

Store in a dry place, away from sunlight.

Store below 25°C.

## **CONTENTS OF THE PACKAGE**

One assembled and pre-filled glass syringe vial, instructions for use, and tracking labels.

The solution contained in the syringe vial is sterile and apyrogenic. The outer surface of the syringe vial is not sterile.

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