EN

MACRO-LOCK POST® ILLUSION® X-RO® FOR DENTAL USE ONLY

DESCRIPTION

MACRO-LOCK POST ILLUSION® X-RO® are dental posts used for reconstruction of devitalized teeth. They are made of X-RO yellow translucent fibers bound in a colored resin Epoxy-vinyl ester matrix, which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray or re-treatment is required.

The **MACRO-LOCK POST ILLUSION® X-RO®** is a serrated-surfaced post, apically tapered with a length of 1.5 mm. The post has fluted slots around its axis for additional macro-retention to the adhesive used, by increasing the surface area available on the post for cement contact.

The radiopacity level of the **MACRO-LOCK POST ILLUSION® X-RO®** is 340% Al eq. (ISO 4049).

The devices are available in 6 sizes and are used in association with an universal pre-shaping drill to perform a first hole and finishing drill associated to each post size according to an O-ring color (See Figure 1). The drills associated to **MACRO-LOCK POST ILLUSION® X-RO®** are **MACRO-LOCK®** drills.

	• 1	• 2	• 3	• 4	• 5	• 6
Ø apical tip	0,80	0,80	1,00	1,00	1,30	1,30
Ø Post head	1,35	1,47	1,67	1,83	2,05	2,20
Starter drill	#1	#1	#2	#2	#2	#2
Finishing Drill	#1	#2	#3	#4	#5	#6

COMPOSITION

MACRO-LOCK POST ILLUSION® X-RO® posts are fiber-reinforced polymer, made of X-RO fiber (80%) and epoxy resin matrix (20%) which contains colored pigment.

MACRO-LOCK® drills are made of stainless steel.

INTENDED USE

MACRO-LOCK POST ILLUSION® X-RO® posts have been designed to support and secure the coronal restoration in case of insufficient residual tooth structure (<4 mm).

MACRO-LOCK® drills are intended to prepare the tooth for insertion of a post that will support and secure the coronal restoration, where insufficient residual tooth structure (<4 mm) exists.

These products should only be used in hospitals, clinics or dental office with qualified dental personnel (dentists).

INDICATIONS FOR USE

MACRO-LOCK POST ILLUSION® X-RO® posts and associated drills are indicated in case of insufficient residual tooth structure (<4 mm) to support and secure the coronal restoration.

CONTRAINDICATIONS

MACRO-LOCK POST ILLUSION® X-RO® posts and associated drills are contraindicated if coronal tooth structure is less than 1.5 mm.

WARNINGS

None known to date.

PRECAUTIONS

There is no data for the use of **MACRO-LOCK POST ILLUSION® X-RO®** posts and associated drills for pregnant women, breastfeeding women and children under 8 years old (for posts only). It is not recommended for such use.

If the device is used in a way different than the one recommended in the IFU, loosening, breakage of the post or extraction of the tooth might occur.

MACRO-LOCK POST ILLUSION® X-RO® posts:

- The post must be immersed or soaked in isopropanol for 1 minute and dried gently with air syringe for 1.5 at least prior to bonding before insertion in the canal. Avoid touching the posts with your fingers after that.

The post can neither be sterilized nor disinfected in thermodesinfector and the use of autoclave is not recommended.

- Posts are single use devices. They must not be re-used to avoid the risk of contamination.

- Shortening the post with a diamond disc should be done outside of the mouth. Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.

- The use of rubber dam is recommended.

- The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

MACRO-LOCK® drills:

- MACRO-LOCK POST ILLUSION® X-RO®** drills are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermodesinfector.

- Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull.

- The universal drill may be used only to remove root canal filling partially, and the finishing drills may be used only to prepare root canal.

- The use of rubber dam is recommended.

- The lifetime of use for the **MACRO-LOCK®** drills stated is 15 uses, following this the dental surgeon is informed that the drill should be replaced.

- Drills must be sterilized individually in an autoclave with the following settings**:

- Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use.

- Autoclave: type B complying with regulation standard in force.

- Sterilization temperature: 134 °C – sterilization time: 18 minutes.

- Store the sterilized components in a dry, dust-free place.

- If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol described in the IFU.

ADVERSE REACTIONS

- Debonding of the post
- Root fracture
- Fracture of the post

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